



SPECIFICATION SHEET

Mayer's Hematoxylin

DESCRIPTION:

Mayer's Hematoxylin is a stain with histological applications. This formula is based on the original formulation developed by Mayer and requires slightly longer application periods.

AVAILABILITY:

<u>CATALOG #</u>	<u>VOLUME</u>
TA-060-MH	60ml
TA-125-MH	125ml

PROTOCOL:

1. Rehydrate tissue if needed.
2. Stain nuclei 1-2 minutes with hematoxylin.
3. Place in Bluing Reagent until stain is blue.
4. Rinse in DI water.
5. Coverslip

RESULTS:

Nuclei are blue.

STORAGE:

Store at room temperature.

PRECAUTIONS:

47790 Westinghouse Dr. Fremont CA 94539 USA
Tel: (800) 828-1628 Fax: (510) 991-2826
Email: labvision@labvision.com Website: www.labvision.com





Do not pipet by mouth. Avoid contact with skin and eyes. Wash after use. Observe all federal, state and local environmental regulations regarding disposal.

LIMITATIONS:

This product is available for research use only; not for diagnostic or therapeutic work. Lab Vision Corporation will not be held responsible for patent infringement or other violation that may occur with the use of this product.

47790 Westinghouse Dr. Fremont CA 94539 USA
Tel: (800) 828-1628 Fax: (510) 991-2826
Email: labvision@labvision.com Website: www.labvision.com

