



SPECIFICATION SHEET

Large Volume Nuclear Fast Red Solution

DESCRIPTION:

Nuclear Fast Red is a stain with histological applications.

AVAILABILITY:

<u>CATALOG #</u>	<u>VOLUME</u>
TA-060-NF	60ml
TA-125-NF	125ml

PROTOCOL:

1. Rinse slides in deionized water.
2. Stain slides for 1 to 5 minutes in Nuclear Fast Red solution.
3. Rinse slides in deionized water.
4. Dehydrate, clear and mount with a non-aqueous mounting medium.

RESULTS:

Nuclei are red, cytoplasm is pale pink.

STORAGE:

Store at room temperature.

PRECAUTIONS:

47790 Westinghouse Dr. Fremont CA 94539 USA
Tel: (800) 828-1628 Fax: (510) 991-2826
Email: labvision@labvision.com Website: www.labvision.com





Do not pipet by mouth. Avoid contact with skin and eyes. Wash after use. Observe all federal, state and local environmental regulations regarding disposal.

LIMITATIONS:

This product is available for research use only; not for diagnostic or therapeutic work. Lab Vision Corporation will not be held responsible for patent infringement or other violation that may occur with the use of this product.

47790 Westinghouse Dr. Fremont CA 94539 USA
Tel: (800) 828-1628 Fax: (510) 991-2826
Email: labvision@labvision.com Website: www.labvision.com

