

Vision Mount

Please note this data sheet has been changed effective December 19, 2011

INTENDED USE

For In Vitro Diagnostic Use

<u>AVAILABILITY:</u>	<u>Catalog #</u>	<u>Volume</u>
	TA-060-UG	60 ml
	TA-125-UG	125 ml

<u>SPECIFICITY:</u>	N/A
<u>ENZYME:</u>	N/A
<u>CHROMOGEN/SUBSTRATE:</u>	N/A

DESCRIPTION

Vision Mount is formulated for covering tissue sections and cell preparations immunohistochemically stained with chromogens producing alcohol soluble OR alcohol insoluble end products. Vision Mount can be used with chromogens such as Fast Red or AEC that require an aqueous mounting medium, or those that do not, such as DAB. No heating is required prior to use. **This product is optimal for use both with and without a glass coverslip.**

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature. Do not use past the expiration date.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

Apply several drops of Vision Mount to stained specimens. Allow to dry.

REFERENCES

N/A

TROUBLESHOOTING

Please contact Thermo Fisher Scientific Technical Support by phone (1-510-991-2800 or 1-800-828-1628) or by email (lab.reagents@thermofisher.com).

