

**Cdk5 Ab-4 (DC17 + DC34)**

Mouse Monoclonal Antibody

**Cat. #MS-1059-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200µg/ml)** (Purified Ab with BSA and Azide)**Cat. #MS-1059-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml)** (Purified Ab without BSA and Azide)**Cat. #MS-1059-PCL (0.1ml)** (Positive Control for Western Blot)**Please note this data sheet has been changed effective December 9, 2011**

**Description:** Cdk5 is expressed in both cycling and noncycling cells and is present in numerous mammalian tissues. Neuronal and muscle cells contain the highest amount of this protein. The p35 protein, which is expressed solely in the brain, activates Cdk5. Cdk5 activity is involved in terminal differentiation of neurons and muscle cells.

**Mol. Wt. of Antigen:** 31kDa

**Epitope:** Not determined

**Species Reactivity:** Human, Mouse and Rat. Others-not known.

**Clone Designation:** DC17 + DC34

**Ig Isotype:** IgG<sub>1</sub> + IgG<sub>1</sub>

**Immunogen:** Recombinant cdk5 protein.

**Applications and Suggested Dilutions:**

- Immunofluorescence
- Immunoprecipitation (Native verified)  
(Use Protein G; Ab 2µg/mg protein lysate)
- Western Blotting (Ab 1-2µg/ml for 2hrs at RT)

The optimal dilution for a specific application should be determined by the investigator.

**Positive Control:** IMR-5 or Ls174T cells.

**Cellular Localization:** Nuclear and cytoplasmic

**Storage and Stability:**

Ab with sodium azide is stable for 24 months when stored at 2-8°C. Antibody WITHOUT sodium azide is stable for 36 months when stored at below 0°C.

**Supplied As:**

200µg/ml antibody purified from the ascites fluid by Protein G chromatography. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide. Also available without BSA and azide at 1mg/ml,

**Suggested References:**

1. Meyerson M, et. al. (1992) EMBO J, 11:2909-2917.
2. Nurse P and Bisset Y. (1981) Nature, 292:558-560.
3. Matsushima H, et. al. (1991) Cell, 65:701-713.
4. Xiong Y, et al. (1992) Cell, 71:505-514.

**Limitations and Warranty:**

Our products are intended FOR RESEARCH USE ONLY and are not approved for clinical diagnosis, drug use or therapeutic procedures. No products are to be construed as a recommendation for use in violation of any patents. We make no representations, warranties or assurances as to the accuracy or completeness of information provided on our data sheets and website. Our warranty is limited to the actual price paid for the product. Lab Vision is not liable for any property damage, personal injury, time or effort or economic loss caused by our products.

**Material Safety Data:**

This product is not licensed or approved for administration to humans or to animals other than the experimental animals. Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. The material contains 0.09% sodium azide as a preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material as indicated above. The National Institute of Occupational Safety and Health has issued a bulletin citing the potential explosion hazard due to the reaction of sodium azide with copper, lead, brass, or solder in the plumbing systems. Sodium azide forms hydrazoic acid in acidic conditions and should be discarded in a large volume of running water to avoid deposits forming in metal drainage pipes.

**For Research Use Only**

