

CD54 / ICAM-1 Ab-1 (Clone 15.2)

Mouse Monoclonal Antibody

Cat. #MS-114-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)**Cat. #MS-114-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml)** (Purified Ab without BSA and Azide)**Please note this data sheet has been changed effective December 6, 2011**

Description: ICAM-1 (intercellular adhesion molecule) has 7 potential N-linked glycosylation sites. It is a single chain glycoprotein of Ig supergene family, present on unstimulated endo-thelial cells (EC) and on a variety of other cell types including activated fibroblasts, EC, macrophages, and lymphocytes. ICAM-1 mediates cell adhesion by binding to integrins CD11a/CD18 (leukocyte adhesion molecule, LFA-1) and to CD11b/CD18 (Mac-1). This interaction enhances antigen-specific T-cell activation. ICAM-1 also binds to CD43 and to *Plasmodium falciparum* infected RBCs.

Comments: Ab-1 inhibits adhesion of infected RBCs and LFA-1 to ICAM-1³.

Mol. Wt. of Antigen: 85-115kDa

Epitope: Not determined

Species Reactivity: Human. Others-not known.

Clone Designation: 15.2

Ig Isotype / Light Chain: IgG₁ / κ

Immunogen: Human monocytes

Applications and Suggested Dilutions:

- Blocks binding of LFA-1 & infected RBCs to ICAM³ (Order Ab without azide)
- Flow Cytometry (For direct, use 10µl of conjugated Ab/10⁶ cells) (For indirect, 0.5µg unconjugated Ab/10⁶ cells)
- Immunohistology (Acetone fixed frozen)

The optimal dilution for a specific application should be determined by the investigator.

Positive Control: Tonsil

Cellular Localization: Cell membrane

Supplied As:

200µg/ml of antibody purified from ascites fluid by Protein G chromatography. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide. Also available without BSA and azide at 1mg/ml.

Storage and Stability:

Ab with sodium azide is stable for 24 months when stored at 2-8°C. Antibody WITHOUT sodium azide is stable for 36 months when stored at below 0°C.

Key References:

1. Bailly P; et al. Proc Natl Acad Sci USA, 1994 7, 91:5306-10.
2. Fujita H; et al. Archives of Biochemistry and Biophysics, 1994, 309:62-9.
3. Berendt AR; et al. Cell, 1992, 68(1):71-81.
4. J Cell Biol (1992) 116:1527-1535.

Limitations and Warranty:

Our products are intended FOR RESEARCH USE ONLY and are not approved for clinical diagnosis, drug use or therapeutic procedures. No products are to be construed as a recommendation for use in violation of any patents. We make no representations, warranties or assurances as to the accuracy or completeness of information provided on our data sheets and website. Our warranty is limited to the actual price paid for the product. Lab Vision is not liable for any property damage, personal injury, time or effort or economic loss caused by our products.

Material Safety Data:

This product is not licensed or approved for administration to humans or to animals other than the experimental animals. Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. The material contains 0.09% sodium azide as a preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material as indicated above. The National Institute of Occupational Safety and Health has issued a bulletin citing the potential explosion hazard due to the reaction of sodium azide with copper, lead, brass, or solder in the plumbing systems. Sodium azide forms hydrazoic acid in acidic conditions and should be discarded in a large volume of running water to avoid deposits forming in metal drainage pipes.

For Research Use Only

CD54 / ICAM-1 Ab-1 (Clone 15.2)

Mouse Monoclonal Antibody

Cat. #MS-114-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)

Cat. #MS-114-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml) (Purified Ab without BSA and Azide)

Please note this data sheet has been changed effective December 6, 2011

Additional Key References:

- Bailly P; Hermand P; Callebaut I; Sonneborn HH; Khamlichi S; Mornon JP; Cartron JP. The LW blood group glycoprotein is homologous to intercellular adhesion molecules. Proceedings of the National Academy of Sciences of the United States of America, 1994 7, 91:5306-10.
- Fujita H; Morita I; Murota S. A possible mechanism for vascular endothelial cell injury elicited by activated leukocytes: a significant involvement of adhesion molecules, CD11/CD18, and ICAM-1. Archives of Biochemistry and Biophysics, 1994, 309:62-9.
- Ganser A; Seipelt G; Verbeek W; Ottmann OG; Maurer A; Kolbe K; Hess U; Elsner S; Reutzel R; Wormann B; et al. Effect of combination therapy with all-trans-retinoic acid and recombinant human granulocyte colony-stimulating factor in patients with myelodysplastic syndromes. Leukemia, 1994, 8:369-75.
- Appleyard RF; Cohn LH. Myocardial stunning and reperfusion injury in cardiac surgery. Journal of Cardiac Surgery, 1993, 8:316-24.
- Bouillon M; Audette M. Transduction of retinoic acid and gamma-interferon signal for intercellular adhesion molecule-1 expression on human tumor cell lines: evidence for the late-acting involvement of protein kinase C inactivation. Cancer Research, 1993, 53:826-32.
- Kojima T; Sasai M; Kobayashi Y. Increased soluble ICAM-1 in tracheal aspirates of infants with bronchopulmonary dysplasia. Lancet, 1993, 342:1023-4.
- Lai CK; Wong KC; Chan CH; Ho SS; Chung SY; Haskard DO; Lai KN. Circulating adhesion molecules in tuberculosis. Clinical and Experimental Immunology, 1993, 94:522-6.
- Liversidge J; Dick A; Cheng YF; Scott GB; Forrester JV. Retinal antigen specific lymphocytes, TCR-gamma delta T cells and CD5+ B cells cultured from the vitreous in acute sympathetic ophthalmitis. Autoimmunity, 1993, 15:257-66.
- Nakayama J; Urabe A; Terao H; Taniguchi S; Hori Y. In situ detection of immunocompetent cells in murine B16 melanoma locally treated with interleukin-2 or microwave hyperthermia. Pigment Cell Research, 1993, 6:111-6.
- Ree HJ; Khan AA; Elsagr M; Liao S; Teplitz C. Intercellular adhesion molecule-1 (ICAM-1) staining of reactive and neoplastic follicles. ICAM-1 expression of neoplastic follicle differs from that of reactive germinal center and is independent of follicular dendritic cells. Cancer, 1993, 71:2817-22.
- Tielemans CL; Delville JP; Husson CP; Madhoun P; Lambrechts AM; Goldman M; Vanherweghem JL. Adhesion molecules and leukocyte common antigen on monocytes and granulocytes during hemodialysis. Clinical Nephrology, 1993, 39:158-65.
- Toyozaki T; Saito T; Takano H; Yorimitsu K; Kobayashi S; Ichikawa H; Takeda K; Inagaki Y. Expression of intercellular adhesion molecule-1 on cardiac myocytes for myocarditis before and during immunosuppressive therapy. American Journal of Cardiology, 1993, 72:441-4.
- Viret C; Davodeau F; Guilloux Y; Bignon JD; Semana G; Breathnach R; Jotreau F. Recognition of shared melanoma antigen by HLA-A2-restricted cytolytic T cell clones derived from human tumor-infiltrating lymphocytes. European Journal of Immunology, 1993, 23:141-6.
- von Stamm U; Brocker EB; von Depka Prondzinski M; Rüter DJ; Rumke P; Broding C; Carrel S; Lejeune FJ. Effects of systemic interferon-alpha (IFN-alpha) on the antigenic phenotype of melanoma metastases. EORTC melanoma group cooperative study No. 18852. Melanoma Research, 1993, 3:173-80.
- Vorkauf M; Duncker G; Nolle B; Sterry W. Adhesion molecules in normal human conjunctiva. An immunohistological study using monoclonal antibodies. Graefes Archive for Clinical and Experimental Ophthalmology, 1993, 231:323-30.
- Yamataka T; Kobayashi H; Yagita H; Okumura K; Tamatani T; Miyasaka M. The effect of anti-ICAM-1 monoclonal antibody treatment on the transplantation of the small bowel in rats. Journal of Pediatric Surgery, 1993, 28:1451-7.
- Fabry Z; Waldschmidt MM; Hendrickson D; Keiner J; Love-Homan L; Takei F; Hart MN. Adhesion molecules on murine brain microvascular endothelial cells: expression and regulation of ICAM-1 and Lgp 55. J of Neuroimmunol, 1992, 36:1-11.
- Morel P; Revillard JP; Nicolas JF; Wijdenes J; Rizova H; Thivolet J. Anti-CD4 monoclonal antibody therapy in severe psoriasis. Journal of Autoimmunity, 1992, 5:465-77.
- Moulton TA; Jiang H; Guarini L; Fetell MR; Fisher PB. Induction of growth suppression and modification of gene expression in multi-drug-resistant human glioblastoma multiforme cells by recombinant human fibroblast and immune interferon. Int J of Cancer, 1992, 51:373-8.



CD54 / ICAM-1 Ab-1 (Clone 15.2)

Mouse Monoclonal Antibody

Cat. #MS-114-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)**Cat. #MS-114-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml)** (Purified Ab without BSA and Azide)

Please note this data sheet has been changed effective December 6, 2011

20. Shimizu Y; Demetris AJ; Gollin SM; Storto PD; Bedford HM; Altarac S; Iwatsuki S; Herberman RB; Whiteside TL. Two new human cholangiocarcinoma cell lines and their cytogenetics and responses to growth factors, hormones, cytokines or immunologic effector cells. *Int J Cancer*, 1992, 52:252-60.
21. Sugama Y; Malik AB. Thrombin receptor 14-amino acid peptide mediates endothelial hyperadhesivity and neutrophil adhesion by P-selectin-dependent mechanism. *Circulation Research*, 1992, 71:1015-9.
22. Sugama Y; Tiruppathi C; offakidevi K; Andersen TT; Fenton JW 2d; Malik AB. Thrombin-induced expression of endothelial P-selectin and intercellular adhesion molecule-1: a mechanism for stabilizing neutrophil adhesion. *Journal of Cell Biology*, 1992, 119:935-44.
23. Taylor PM; Rose ML; Yacoub MH; Pigott R. Induction of vascular adhesion molecules during rejection of human cardiac allografts. *Transplantation*, 1992, 54:451-7.

