

Grb2 / sem5 / ASH Signalling Protein Ab-4 (Clone 2GB04; same C1.4)**Mouse Monoclonal Antibody**

Cat. #MS-1457-P0, -P1, or -P (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)

Cat. #MS-1457-P1ABX or -PABX (0.1ml or 0.2ml at 1.0mg/ml) (Purified Ab without BSA and Azide)

Cat. #MS-1457-R7 (7.0ml) (Ready-to-Use for Immunohistochemical Staining)

Cat. #MS-1457-PCS (5 Slides) (Positive Control for Histology)

Description: Grb2 protein contains one SH2 domain flanked on both sides by an SH3 domain. Ab-1 binds to the junction between the amino terminal SH3 domain and the SH2 domain of Grb2. Reactivity with Grb2 deletion mutants suggests that its epitope maps between aa 54-62. In response to EGF, Grb2 binds, through its SH3 domains, to the guanine nucleotide exchange factor, Sos. The Grb2-Sos complex then binds the EGF-receptor. These interactions play a critical role in the activation of the p21^{ras} signalling pathway.

Mol. Wt. of Antigen: 24-25kDa

Epitope: aa 54-62

Species Reactivity: Human, Cow, Mouse, Rat, and *Xenopus*. Others-not known.

Clone Designation: 2GB04 (same as C1.4)

Ig Isotype / Light Chain: IgG₁ / κ

Immunogen: Purified recombinant human Grb2 protein.

Applications and Suggested Dilutions:

- Immunofluorescence
- Immunohistology (Formalin/paraffin)
(Use Ab at 1-2µg/ml for 30 min at RT)
- * [Staining of formalin-fixed tissues REQUIRES boiling tissue sections in 10mM citrate buffer, pH 6.0, (**NEOMARKERS'** Cat. #AP-9003), for 10-20 min followed by cooling at RT for 20 min.]

The optimal dilution for a specific application should be determined by the investigator.

Positive Control: HeLa cells. Tonsil.

Cellular Localization: Cytoplasmic

Storage and Stability: Ab with sodium azide is stable for 24 months when stored at 2-8°C. Antibody WITHOUT sodium azide is stable for 36 months when stored at below 0°C.

Supplied As:

200µg/ml antibody purified from the ascites fluid by Protein G chromatography. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide. Also available without BSA and azide at 1mg/ml.

or

Prediluted antibody which is ready-to-use for staining of formalin-fixed, paraffin-embedded tissues.

Suggested References:

1. Boulikas T. . Critical Reviews in Eukaryotic Gene Expression, 1995, 5(1):1-77.
2. Chardin P; et al. Febs Letters, 1995 Aug 1, 369(1):47-51.
3. Downward J. Febs Letters, 1994 Jan 31, 338(2):113-7.
4. Schlessinger J. Current Opinion in Genetics and Development, 1994 Feb, 4(1):25-30.

Limitations and Warranty:

Our products are intended FOR RESEARCH USE ONLY and are not approved for clinical diagnosis, drug use or therapeutic procedures. No products are to be construed as a recommendation for use in violation of any patents. We make no representations, warranties or assurances as to the accuracy or completeness of information provided on our data sheets and website. Our warranty is limited to the actual price paid for the product. NeoMarkers is not liable for any property damage, personal injury, time or effort or economic loss caused by our products.

Material Safety Data:

This product is not licensed or approved for administration to humans or to animals other than the experimental animals. Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. The material contains 0.09% sodium azide as a preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material as indicated above. The National Institute of Occupational Safety and Health has issued a bulletin citing the potential explosion hazard due to the reaction of sodium azide with copper, lead, brass, or solder in the plumbing systems. Sodium azide forms hydrazoic acid in acidic conditions and should be discarded in a large volume of running water to avoid deposits forming in metal drainage pipes.

For Research Use Only