

SPECIFICATION SHEET

BCIP/NBT

DESCRIPTION:

BCIP/NBT (5-Bromo-4-chloro-3-Indolyl Phosphate/Nitroblue Tetrazolium) is a refrigerator-stable, single component solution for immunohistochemical, *in situ* hybridization, Western, Northern and Southern blot procedures. This reagent is specifically formulated for use with Alkaline Phosphatase enzyme systems. BCIP/NBT produces a purple reaction product which is soluble in alcohol.

AVAILABILITY:

<u>CATALOG #</u>	<u>VOLUME</u>
TA-060-AB	60ml
TA-125-AB	125ml

IMMUNOHISTOCHEMICAL STAINING PROTOCOL:

1. Apply solution to tissue section after incubation with Alkaline Phosphatase label.
2. Incubate tissue section for 5-15 minutes at room temperature, depending on the desired stain intensity.
3. Rinse tissue, counterstain, and coverslip.

STORAGE:


Store at 2-8°C. Product is stable for 18 months from date of manufacture.

PRECAUTIONS:

47790 Westinghouse Dr. Fremont CA 94539 USA
Tel: (800) 828-1628 Fax: (510) 991-2826
Email: labvision@labvision.com Website: www.labvision.com



TA-060-AB
Rev 21127
Page 1 Of 1



Do not pipet by mouth. Avoid contact with skin and eyes. Wash after use. Observe all federal, state and local environmental regulations regarding disposal.

LIMITATIONS:

This product is available for research use only; not for diagnostic or therapeutic work. Lab Vision Corporation will not be held responsible for patent infringement or other violation that may occur with the use of this product.